


## HoloMonitor® App Suite

# Cell Quality Control

It is well known that when cells are cultured over longer or shorter time periods they may undergo morphological changes. Quick assessment of the quality of a cell culture prior to experiments may therefore be useful to avoid unexpected or inferior results.

### DESCRIPTION

The HoloMonitor® Cell QC Assay presents end-point data on a cell population level in terms of basic morphology, cell count and cell confluence. The assay can be used to ensure that the cells are viable and cell count and/or confluence is appropriate at the start of an experiment. It can also be used as a tool to detect undesired morphological changes in your cell culture when compared to previous experiments.



| Well | Cell count [#/cm <sup>2</sup> ]             | Confluence [%] | Mean cell volum...         | Mean cell area [μ...       |
|------|---|----------------|----------------------------|----------------------------|
| A1   | 3,3×10 <sup>4</sup> (±0)                    | 16 (±0)        | 1,6×10 <sup>3</sup> (±0)   | 4,4×10 <sup>2</sup> (±0)   |
| A2   | 3,3×10 <sup>4</sup> (±3,6×10 <sup>2</sup> ) | 16 (±0,2)      | 1,6×10 <sup>3</sup> (±6,6) | 4,3×10 <sup>2</sup> (±6,8) |
| B1   | 3,3×10 <sup>4</sup> (±0)                    | 16 (±0)        | 1,6×10 <sup>3</sup> (±0)   | 4,2×10 <sup>2</sup> (±0)   |
| B2   | 3,3×10 <sup>4</sup> (±0)                    | 16 (±0)        | 1,6×10 <sup>3</sup> (±0)   | 4,2×10 <sup>2</sup> (±0)   |
| C1   | 3,3×10 <sup>4</sup> (±0)                    | 16 (±0)        | 1,6×10 <sup>3</sup> (±0)   | 4,2×10 <sup>2</sup> (±0)   |
| C2   | 3,3×10 <sup>4</sup> (±0)                    | 16 (±0)        | 1,6×10 <sup>3</sup> (±0)   | 4,2×10 <sup>2</sup> (±0)   |
| D1   | 3,3×10 <sup>4</sup> (±0)                    | 16 (±0)        | 1,6×10 <sup>3</sup> (±0)   | 4,2×10 <sup>2</sup> (±0)   |
| D2   | 3,3×10 <sup>4</sup> (±0)                    | 16 (±0)        | 1,6×10 <sup>3</sup> (±0)   | 4,2×10 <sup>2</sup> (±0)   |
| A3   | 3,3×10 <sup>4</sup> (±0)                    | 16 (±0)        | 1,6×10 <sup>3</sup> (±0)   | 4,2×10 <sup>2</sup> (±0)   |
| A4   | 3,3×10 <sup>4</sup> (±0)                    | 16 (±0)        | 1,6×10 <sup>3</sup> (±0)   | 4,2×10 <sup>2</sup> (±0)   |
| B3   | 3,3×10 <sup>4</sup> (±0)                    | 16 (±0)        | 1,6×10 <sup>3</sup> (±0)   | 4,2×10 <sup>2</sup> (±0)   |
| B4   | 3,3×10 <sup>4</sup> (±0)                    | 16 (±0)        | 1,6×10 <sup>3</sup> (±0)   | 4,2×10 <sup>2</sup> (±0)   |
| C3   | 3,3×10 <sup>4</sup> (±0)                    | 16 (±0)        | 1,6×10 <sup>3</sup> (±0)   | 4,2×10 <sup>2</sup> (±0)   |
| C4   | 3,3×10 <sup>4</sup> (±0)                    | 16 (±0)        | 1,6×10 <sup>3</sup> (±0)   | 4,2×10 <sup>2</sup> (±0)   |
| D3   | 3,3×10 <sup>4</sup> (±0)                    | 16 (±0)        | 1,6×10 <sup>3</sup> (±0)   | 4,2×10 <sup>2</sup> (±0)   |
| D4   | 3,3×10 <sup>4</sup> (±0)                    | 16 (±0)        | 1,6×10 <sup>3</sup> (±0)   | 4,4×10 <sup>2</sup> (±0)   |
| A5   | 3,3×10 <sup>4</sup> (±0)                    | 16 (±0)        | 1,6×10 <sup>3</sup> (±0)   | 4,2×10 <sup>2</sup> (±0)   |
| A6   | 3,3×10 <sup>4</sup> (±0)                    | 16 (±0)        | 1,6×10 <sup>3</sup> (±0)   | 4,2×10 <sup>2</sup> (±0)   |
| B5   | 3,3×10 <sup>4</sup> (±0)                    | 16 (±0)        | 1,6×10 <sup>3</sup> (±0)   | 4,2×10 <sup>2</sup> (±0)   |

### Output - End-point data given as values/well

- Cell area (μm<sup>2</sup>)
- Mean cell volume (μm<sup>3</sup>)
- Cell count (no of cells/well)
- Cell confluence (% cell covered area)

### HOLOMONITOR APP SUITE

HoloMonitor® App Suite is a completely new proprietary software for analysis of images and data generated by the HoloMonitor® M4 base unit. HoloMonitor® App Suite focuses on biological applications and enables researchers within all levels of cell biology to easily perform live-cell studies on various cellular events.

### FURTHER INFORMATION

phiab.se  
info@phiab.se