

# HoloMonitor<sup>®</sup> M4

## CELL MOTILITY PROTOCOL

### MATERIAL


- **HoloMonitor<sup>®</sup> M4**, placed inside a cell incubator.
- **Hstudio software**, version 2.7.1 or later.
- **Culture vessel** by choice:
  - Sarstedt TC Dish 35, Standard (cat. # 83.3900)
  - Sarstedt TC Plate 6 Well, Standard, F (cat. # 83.3920.005)
  - Sarstedt lumox<sup>®</sup> multiwell, 24 Well (cat. # 94.6000.014)
  - Sarstedt lumox<sup>®</sup> multiwell, 96 Well (cat. # 83.3924.005)
- PHI **HoloLid** for selected vessel:

<u>Vessel</u>	<u>HoloLid</u>
Sarstedt TC Dish 35, Standard	71110
Sarstedt TC Plate 6 Well, Standard, F	71120
Sarstedt lumox <sup>®</sup> multiwell, 24 Well	71130
Sarstedt lumox <sup>®</sup> multiwell, 96 Well	71140

HoloLid product information and protocol is available [here](#).

- PHI **Vessel holder** for the selected vessel For information regarding vessel holders contact PHI at [support@phiab.se](mailto:support@phiab.se).
- **Cells** suspended to reach a confluence of 2-5 % when seeded (approx. 6 000-11 000 cell/cm<sup>2</sup> for L929, A375, and Jimt-1 cells). Other cell types may require a different seeding concentration.
- [Setup and Operation Manual](#) for using HoloMonitor M4, if the user is unfamiliar with the imaging procedures.

We recommend that each experiment is performed at least three times to acquire sound statistics.



## PREPARATION

1. Seed the cells to a confluence of 2-5 % in preferred vessel. The final working volumes, essential for using **HoloLids** are:

Volume	Vessel
3.0 ml	Sarstedt TC Dish 35, Standard
3.0 ml/well	Sarstedt TC Plate 6 Well, Standard, F
1.8 ml/well	Sarstedt lumox® multiwell, 24 Well
170 µl/well	Sarstedt lumox® multiwell, 96 Well


2. Put the vessel into the cell incubator and let the cells attach for 2-24h.
3. Sterilize the **HoloLids** according to the **HoloLid protocol**.
4. Add the treatment, if stated in the experimental setup.
5. Put on the standard lid.

## IMAGING

Start up the **HoloMonitor** and proceed with the calibration. The values achieved should lie within the green area of the calibration results bar.


### For imaging with a motorized stage

1. Wipe off the **Vessel holder** with alcohol and put it in to the LAF-bench, the grips facing down.
2. Place the cell samples on to the **Vessel holder**:



3. Replace the standard lids with the appropriate **HoloLid**, following the HoloLid protocol.

- Thereafter place the **Vessel holder** with the samples on the HoloMonitor stage:



- Go to the **Live capture tab** in the Hstudio software and select the appropriate vessel template:

<u>Vessel</u>	<u>Template</u>
Sarstedt TC Dish 35, Standard	Petri dishes 40 mm
Sarstedt TC Plate 6 Well, Standard, F	Sarstedt 6 well plate with PHI lid
Sarstedt lumox® multiwell, 24 Well	Lumox 24 well plate
Sarstedt lumox® multiwell, 96 Well	Lumox 96 well plate with PHI lid

- Create a **Project** for image storage.
- Focus the images at a position close to the center of the plate/vessel.
- Check **Timelapse** and type the total time and interval of the time-lapse imaging. 5 minutes between captures is recommended. However, for very fast cells an even shorter interval may be required.
- Check **Capture pattern** and select the wells and positions to be captured, as described in the [Setup and Operation Manual](#). Alternatively select capture positions and click the **Remember button** for as many positions as required. However, short interval between image captures (< 5 min) limits the number of positions that can be captured.
- Click **Advanced setup** and check **Multiple destination groups**.
- Make sure **One group per well** is not checked.
- Click **Save and close**.
- Click **Capture**.
- Go to the **View image tab** and review the images for quality.
- Await the time-lapse capturing to finish.

## For imaging with a fixed stage


1. Replace the standard lid with the appropriate **HoloLid**, following the HoloLid protocol.
2. Place the sample on the **HoloMonitor stage** using the appropriate distance plate:

<u>Vessel</u>	<u>Distance plate</u>
Sarstedt TC Dish 35, Standard	# 2
Sarstedt TC Plate 6 Well, Standard, F	None
Sarstedt lumox® multiwell, 24 Well	# 1
Sarstedt lumox® multiwell, 96 Well	None

3. Go to the **Live capture tab** and ensure that the images are well focused. Adjust the software focus, if required.
4. Create a **Group**.
5. Check **Timelapse** and type the total time and interval of the time-lapse capture. 5 minutes between captures is recommended. However, for very fast cells an even shorter interval may be required.
6. Click **Capture**.
7. Go to **View image** and review the images for quality.
8. Await the time-lapse capturing to finish.

## ANALYSIS

1. Go to the **Identify cells tab**. Check the segmentation and, if needed, adjust the **Threshold** and **Min object size** settings to fit the cells for all **Groups**. The settings can be applied for all images within each **Group** but need to be validated and possibly adjusted for all **Frames**. Discard bad frames.
2. Go to **Track cells**. Add all frames from one position to the tracking analysis. Individual cells to be tracked are added by clicking on them.
3. Move the **Timeline slider** to the right to see the tracks of the added cells. Adjust the possible errors of the software using the **Warnings list**.
4. When the cells divide, a warning sign appears in the frame. This indicates a cell division. Click **Division** if both daughter cells are to be tracked (optional). Otherwise, a random daughter cell is selected and subsequently tracked by the software.
5. Go to **Plot movement** and check all cells to be included in the plot. The colored tracks show the movements of the selected cells, with the origin as the starting point for each cell:



The plot can be saved as an image in several formats. If different samples (treatments) are to be compared it is recommended to adjust the X- and Y-scales to be identical for all samples before export to XML-files.

