

HOLOMONITOR® APP SUITE PROTOCOL

WOUND HEALING ASSAY

This protocol helps to set up a **Wound Healing Assay** using HoloMonitor® M4 and the HoloMonitor® App Suite software. The HoloMonitor® Wound Healing Assay provides an automated label-free wound healing assay, measuring gap closure. Additionally, selected cells can be individual cell tracked for detailed individual cell movement and morphology analysis.

REQUIREMENTS:

- HoloMonitor® M4, placed in incubator
- HoloMonitor® M4 App Suite
- Culture vessel of choice with cells
- HoloLid™ for selected vessel
- Vessel holder for selected vessel


OUTPUT:

- Gap width (µm)
- Cell covered area (µm² and %)
- Cell free area (µm² and %)

REANALYSIS:


- **Guided assays**
 - Cell QC
 - Cell Proliferation
 - Cell Motility
 - Dose Response
- **In-depth assays**
 - Spatial tracking
 - Wound Healing
 - Cell Morphology

WORK FLOW:






PREPARATIONS

Materials


- ✓ HoloMonitor® M4, placed inside the incubator
- ✓ HoloMonitor® App Suite software
- ✓ Cell culture vessel.
 - ▶ We recommend using ibidi® μ-Dish 35 mm with Culture-Insert 2 Well high (cat. # 81176) or ibidi® μ-Plate 24-well with Culture-Insert 2 Well (cat. # 80241)
- ✓ HoloLid™ for the selected vessel
- ✓ Vessel holder for the selected vessel
- ✓ Cells
- ✓  Setup and Operational Manual for HoloMonitor® M4

Steps

1. Seed the cells to at least **90 % confluence** in the inserts according to  ibidi protocol.
2. Place the vessel into the incubator and let cells attach for 2-24 hours.
3. **Start the software** and wait for complete instrument initialization.
4. Run an **auto-calibration**. With successful calibration, the instrument is ready to use.
5. **Sterilize the HoloLids™** according to the specific  HoloLid™ protocol.
6. Remove the inserts and add respective cell **treatment**. The **final working volume** per well for recommended ibidi® vessels is 2.5 mL.
7. **Slide the cell culture vessel** onto the Vessel holder, its grips facing towards you. Ensure that the vessel is parallel to the holder. There is a spring that holds the vessel in place.
 - ▶ When using multi-well plates, place them with the cut-off corner to the left.
8. **Replace the standard lids with the HoloLid™**.
9. **Put the vessel holder with the sample on the Holo-Monitor® M4 stage** and click it to secure.
10. Select the **Wound Healing Assay** and proceed by clicking the **Setup Assay** button.



Successful auto-calibration window




AppSuite main window with selected Wound Healing Assay

EXPERIMENT SETUP


1 Basic setup: describe the experiment and assign treatments/conditions to the wells


1. Enter the experiment **name**, optional experiment **description** and cell types.
2. Select the correct **vessel map** from the drop-down list.
3. Map **treatments and conditions** on the vessel map. **Select wells** by marking them with the left mouse button while moving the cursor over the relevant well/s.
4. Add the **treatment name/s** in the text box below the vessel map and click **Add** /press Enter. It is possible to add wells as **individual treatments**. Marked well/s are light blue, selected wells will appear dark blue.
5. Proceed to **Capture setup**.




Basic Setup window

2 Capture setup: Select the experiment time settings, choose capture positions

1. Adjust the default settings for **duration**, **interval**.
2. Add **capture positions**: The position list is open by default. Click positions on the vessel map and add them to the position list with the **Add current location button**. In case the image quality is poor, a warning sign  appears. **Adjust focus or position location** if necessary.
▶ Note that the gap might not be exactly where the vessel map indicates.
3. Ensure that the **storage requirement** for the experiment does not exceed the computer capacity.
4. Run a full or quick **validation** of the selected positions to ensure **good image quality**.
5. When satisfied with the experiment setup, click **Proceed to Capture**.




Capture Setup window

MANUAL FOCUS ADJUSTMENT: The focusing tool is located in the **Controls** tab. Move the black square or click on the Arrow buttons to move the stage up and down. Save an adjusted focus setting for the selected position by using the **Apply Current Focus** button. For details, consult the  Setup and Operational Manual.

3 Capture: Review the experiment in real-time during the time-lapse

1. Click **Start Capture**.
2. To stop the experiment ahead of time, click the stop button. Note that it is **NOT** possible to restart the experiment once it has been stopped.
3. Go to the **Experiments** tab and open your ongoing experiment to preview the captured images during the run.
▶ Wait for the experiment to finish before starting In-depth Analysis.
4. When the Experiment capture finished, click the **Show Result** button to get directly to the **Results** page.




Capture window

RESULTS & ANALYSIS

Experiments tab


1. Click Experiments to see a list of the experiments.
2. Click on the experiment title to open an experiment summary.
3. Click Open to open the results page.



Experiments tab

Experiment overview tab


1. See the experiment summary, view all images and go to the experiment setup by choosing the respective tab.
2. Generate in-depth analysis data from the captured images by clicking on the Wound Healing icon. A new window for the in-depth analysis will open.
3. Create New Guided Assay Results from this experiment by clicking the respective button.




Experiment overview tab

In-depth analysis — wound healing tab


1. Begin by adding frames to the analysis either by drag-and-drop or using the Add selected or Add all button.
 - ▶ Check the image quality before using it for in-depth analysis. See the [image guide](#) for more information.
2. Adjust the image threshold and object size and click Apply to all frames. You can alter the viewing options, too.
3. Based on the identification in the previous step, result values are generated and displayed in the table:
 - ▶ gap with (mm), cell-covered area (% and μm^2) and cell-free area (% and μm^2).
4. Check the gap width values in the table to ensure that the values decrease evenly with time.
 - ▶ In case of outliers, select the outlier frame in the list and adjust the threshold for that specific frame.
5. Export to Excel for further analysis. The exported data include: gap width (μm), cell-covered area (% and μm^2) and cell-free area (% and μm^2) and the settings overview.



The HoloMonitor® App Suite experiment data can be re-analyzed to generate results for other assay types. Find more information in the [overview](#) and the [Setup and Operational Manual](#).



Wound healing tab



Gap closure graph

Calculate the cell front velocity by dividing the first gap width value of the linear phase of the graph by the time value at the end of the slope's linear phase. In the example above: $410 \mu\text{m} / 10.5 \text{ h} = 39 \mu\text{m/h}$